





# 20 PRODUCT 21 CATALOG

America's Dental  
Equipment Company



(800) 423-5657 VISIT US AT [WWW.DHPDENTAL.COM](http://WWW.DHPDENTAL.COM)





# America's Dental Equipment Company

## For Over 50 Years

### Dansereau Health Products, Inc.

Dansereau Health Products, Inc has been manufacturing and selling direct to the dentist quality dental equipment for over 50 years. All our manufacturing is done in our factory in Corona, CA USA. Many of our competitors will import parts and do a modest assembly and consider their product manufactured in the USA.

Dansereau has a complete machine shop that fabricates all our parts in our factory in Corona, CA USA. We use our fabricated parts in the assembly of all our products, allowing us the most stringent of quality controls resulting in a 10 year warranty, the best in the dental industry.

### Warranty

Dansereau includes with all their equipment one of the most comprehensive and complete warranties. 10 Years on our Dental Chairs and Delivery Systems.

### Support

Dansereau has designed a complete system to support your dental equipment for decades to come. We have approximately 300 independent service technicians available for field support. We have designed a series of video's to instruct on simple maintenance of your dental equipment and have a comprehensive parts website to easily search for any maintenance and service parts.

### Technology Support

Dansereau has moved Customer Service into the 21st Century with our innovative approach to supporting the dentist in their use of Dansereau Dental Equipment. Video conferencing, Online Video are a few of Dansereau's methods of supporting dentists, office staff and service technicians.






**Industry Leading 10 Year Warranty**  
On All Dental Chairs



Over 150 Color Choices  
Call for color samples



Memory Foam Cushions



Dual Hand's Controls  
Operation Auto Return  
& Adjustable Preset



Drop Down Arms for Easy Access



800 lb Lift Capacity



**\$2995**  
- Ultra Thin Seat Back  
- Aseptic Cushions  
- 150 Color Choices  
- 10 Year Warranty

AM3 Traditional  
Dental Chair



**\$3395**  
- Ultra Thin Seat Back  
- Aesthetic Memory Foam  
Cushions  
- 150 Color Choices  
- 10 Year Warranty

AM3 Contemporary  
Dental Chair



**\$3495**  
- Ultra Thin Seat Back  
- Memory Foam Cushions  
- Increased Chair Height 35"  
- Sedation Surgery Arm Supports  
- 150 Color Choices  
- 10 Year Warranty

AM Surgery  
Dental Chair



**\$2795**  
- Ultra Thin Seat Back  
- Removable/Reversible Cushions  
- 150 Color Choices  
- 10 Year Warranty

California  
Dental Chair



**\$2995**  
- Ultra Thin Seat Back  
- Pediatric Seat Cushion  
- Narrower / Shorter  
- 150 Color Choices  
- 10 Year Warranty

AM Pedo  
Dental Chair



**\$3295**  
- Ultra Thin Seat Back  
- Memory Foam Cushions  
- Extreme Patient Comfort  
- 150 Color Choices  
- 10 Year Warranty

Infinity  
Dental Chair





## AMERICA CART OPERATOR 10 Year Warranty

- AM3 Traditional Dental Chair
- 3 HPC Unit and 3 Way Syringe
- U Cart with Adjustable Height
- Complete Assistants Package
- Diamond Dental Light
- 800 lb Lift Capacity
- Over 150 Color Choices

\$5395



## UNIVERSAL CART OPERATOR 10 Year Warranty



- AM3 Traditional Dental Chair
- Radius 3 HPC Unit and 3 Way Syringe
- Radius Telescoping Assistants Package
- U-Cart with Adjustable Height
- Laminated Worksurface Area
- Diamond Dental Light
- 800 lb Lift Capacity
- Over 150 Color Choices

\$5595







## AMERICA OPERATOR 10 Year Warranty

- AM3 Traditional Dental Chair
- 3 HPC Unit and 3 Way Syringe
- Telescoping Assistants Package
- Diamond Dental Light
- 800 lb Lift Capacity
- Over 150 Color Choices

**\$5895**

Includes Free Shipping\*



## UNIVERSAL OPERATOR 10 Year Warranty



- AM3 Traditional Dental Chair
- Radius 3 HPC Unit and 3 Way Syringe
- Radius Telescoping Assistants Package
- Diamond Dental Light
- 800 lb Lift Capacity
- Over 150 Color Choices

**\$6795**



Terms Apply\* 48 Continental USA  
- Out of Area Shipments May Incur Modest Additional Costs





## INDEPENDENCE OPERATORY 10 Year Warranty

- Infinity Dental Chair
- 3 HPC Unit and 3 Way Syringe
- 3 Way Syringe w/ Assistants Package
- Height Adjustable Work Surface
- Radius Diamond Dental Light
- 800 lb Lift Capacity
- Over 150 Color Choices

\$6195



## ORTHODONTIC DENTAL CHAIRS Quality Materials, Quality Products



\$2995

### Classic Ortho Chair

- Power Base, Powered Seat Back
- Lift Up Arms, Synchronized Arm Movement,
- Thin Seat Back, Foot Operated,
- 150 Color Choices
- 10 Year Warranty



\$1995



### Simple Ortho Chair

- Fixed Base, Powered Seat Back,
- Lift Up Arms, Synchronized Arm Movement,
- Thin Seat Back, Foot Operated,
- 150 Color Choices
- 10 Year Warranty

Monitor Bracket  
\$195

Minitrol W/ Asst  
Package \$1495

120V Convenience  
Outlet \$125



Orthodontic Cabinet  
\$1295



12 O'CLOCK DELIVERY SYSTEMS



INDEPENDENCE  
CABINET/WALL SYSTEM

\$2495



UNIVERSAL CART SYSTEM

\$2495



INFINITY CART SYSTEM

\$1895

- 3 HPC Automatic Unit w/ 3-Way Syringe
- Complete Assistants Package w/ HVE and Saliva Ejector
- Height Adjustable Work-surface
- Wall or Cabinet Installation
- Complete Junction Box (Not Shown)
- 10 YEAR WARRANTY

\*Water Bottle System Not Included

- 3 HPC Automatic Unit w/ 3-Way Syringe
- Complete Assistants Package w/ HVE and Saliva Ejector
- Height Adjustable U Cart w/ work surface and EZ Roll Casters
- Unit Foot Control
- Complete Junction Box (Not Shown)
- 10 YEAR WARRANTY

- 3 HPC Automatic Unit w/ 3-Way Syringe
- Complete Assistants Package w/ HVE and Saliva Ejector
- Height Adjustable U Cart w/ EZ Roll Casters
- Unit Foot Control
- Flex Arm with Air Lock
- Complete Junction Box (Not Shown)
- 10 YEAR WARRANTY

\*Water Bottle System Not Included

CABINET MOUNTED DELIVERY SYSTEMS

INFINITY FLEX UNIT | \$2495



- 3 Handpiece Automatic System
- Dentists 3-Way Syringe
- Straight Soft Tubing
- Assistants Package
- High Evacuation Valve
- Saliva Ejector
- Assistants 3-Way Syringe
- Complete Dental Junction Box

\*Water Bottle System Not Included

INFINITY FACEPLATE UNIT | \$1795




- 2 Handpiece Automatic System
- Dentists 3-Way Syringe
- Straight Soft Tubing
- Assistants Package
- High Evacuation Valve
- Saliva Ejector
- Assistants 3-Way Syringe
- Complete Dental Junction Box

INFINITY SLIDE OUT UNIT

\$1895 (Without Asst. Package)

\$2195 (With Asst. Package)



- 3 Handpiece Automatic System
- Dentists 3-Way Syringe
- Straight Soft Tubing
- Cabinet Mounted Asst. Package
- High Evacuation Valve
- Saliva Ejector
- Assistants 3-Way Syringe
- Complete Dental Junction Box

\*Water Bottle System Not Included





OPERATOR SEATING



SIGNATURE ASSISTANT'S  
STOOL

- Adjustable Foot Ring
- Aseptic Cushions
- Aluminum Star
- Chrome Cylinder
- Non - Marr Casters

\$625



SIGNATURE DENTIST'S  
STOOL

- Ergonomic Controls
- Aseptic Cushions
- Aluminum Star
- Chrome Cylinder
- Non - Marr Casters

\$475



SIGNATURE SADDLE  
ASSISTANT'S STOOL

- Adjustable Foot Ring
- Aseptic Cushions
- Aluminum Star
- Chrome Cylinder
- Non - Marr Casters

\$725



CALIFORNIA DENTIST'S  
STOOL

- Single Lever Lift
- Aseptic Cushions
- Poly Resin Star
- Black Cylinder
- Non - Marr Casters

\$395



BASIC DENTIST'S  
STOOL

- Single Lever Lift
- Aseptic Cushions
- Aluminum Star
- Chrome Cylinder
- Non - Marr Casters

\$275



SADDLE DENTIST'S  
STOOL

- Aseptic Cushions
- Aluminum Star
- Black Cylinder
- Non - Marr Casters

\$395


OPERATOR SEATING



INFINITY ASSISTANTS SYSTEMS

Installs on Any Chair

INFINITY DELIVERY SYSTEMS



INDEPENDENCE  
ASSISTANTS SYSTEM

- Radius Assistants Package
- Height Adjustable Laminated Table
- Assistants 3-Way Syringe
- 2 HVE Valves with 1 Saliva Ejector
- Water Quick Disconnect
- 10 YEAR WARRANTY

\$1495

INFINITY 1  
ASSISTANTS PACKAGE

- Full Telescoping Assistants Arm
- Post Mounting Bracket
- Assistants 3-Way Syringe
- 2 HVE Valves with 1 Saliva Ejector
- Water Quick Disconnect
- 10 YEAR WARRANTY

\$995

INFINITY 2  
TELESCOPING ARM  
PACKAGE

- Short Telescoping Assistants Arm
- Post Mounting Bracket
- Assistants 3-Way Syringe
- 2 HVE Valves with 1 Saliva Ejector
- Water Quick Disconnect
- 10 YEAR WARRANTY

\$695

CUSPIDOR  
SYSTEM

- Cuspidor - Ceramic Bowl
- Short Telescoping Assistants Arm
- Timed Rinse Feature
- Manual Cupfiller
- Assistants 3-Way Syringe
- 2 HVE Valves with 1 Saliva Ejector
- Water Quick Disconnect
- 10 YEAR WARRANTY

\$1595

INFINITY 1  
DELIVERY SYSTEM

- 3 Handpiece Automatic System
- Dentists 3-Way Syringe Included
- Recessed Tray
- Straight Soft Tubing
- Flex Arm with Air Lock
- Telescoping Assistants Arm
- Assistants 3-Way Syringe
- 2 HVE Valves with Saliva Ejector
- Water Quick Disconnect
- Mounting Post and Complete Junction Box

\$2595

INFINITY 2  
DELIVERY SYSTEM

- 3 Handpiece Automatic System
- 3-Way Syringe Included
- Recessed Tray
- Straight Soft Tubing
- Flex Arm with Air Lock
- 10 YEAR WARRANTY
- Unit Mounted, Cabinet Mounted or Wall Mounted Available

\$1795

INFINITY 3  
DELIVERY SYSTEM

- 3 Handpiece Automatic System
- Dentists 3-Way Syringe Included
- Recessed Tray
- Straight Soft Tubing
- Flex Arm with Air Lock
- Cuspidor with Timed Rinse Bowl
- Manual Cupfiller
- Telescoping Assistants Arm
- Assistants 3-Way Syringe
- 2 HVE Valves with Saliva Ejector
- Water Quick Disconnect
- Mounting Post and Complete Junction Box

\$3195



# AIR COMPRESSOR SYSTEMS



## QUIET FLOW 1HP AIR COMPRESSOR

- 1.25 HP Oil Less Motor
- 120 Volt or 240 Volt
- Regenerating Drying System
- Automatic Purge
- 99.9% Coalescing Filter
- Installation Hardware Included
- Quiet Operation 40db
- Limited 5 Year Warranty

\$3795



## QUIET FLOW 2 DUAL HP AIR COMPRESSOR

- 2.5 HP Oil Less Motor
- 240 Volt
- Regenerating Drying System
- Automatic Purge
- 99.9% Coalescing Filter
- Installation Hardware Included
- Quiet Operation 40db
- Limited 5 Year Warranty


\$4895



CLEAN FLOW 3.75 HP  
(TRIPLE HEAD)  
COMPRESSOR WITH DRYER  
\$6795



CLEAN FLOW 5 HP  
(QUAD HEAD)  
COMPRESSOR WITH DRYER  
\$7995



# COMPLETE DENTAL OFFICES



2 Complete Operatories With  
Stools & X-Ray Machine



Autoclave Sterilizer

## Complete Utility Room

- 1 HP Oil Less
- Compressor with Dryer
- 1 HP Vacuum System
- Air Water Separator
- Low Voltage Control Box
- Low Voltage Control Panel

\$26,995



Consider a Complete Dental Office Package when you are planning your new Dental Office.

All New Equipment, Fully Warrantied and Manufactured in the USA.

Many other packages are available, ask your Dansereau Rep for more information.



# WET VACUUM SYSTEMS



## \$1795 1HP PRECISION VACUUM SYSTEM

- 1 HP Solid State Motor
- 120Volt or 240Volt (Specify on Order)
- Machined Naval Brass Housing
- Spin Balance Impeller for Longevity
- Particulate Filter Included
- Installation Hardware Included
- Add an Air/Water Separator \$199
- Limited 5 Year Warranty



## \$2395 NEW! 2HP PRECISION VACUUM SYSTEM

- 2 HP Solid State Motor
- 240 Volt (Specify on Order)
- Machined Naval Brass Housing
- Spin Balance Impeller for Longevity
- Particulate Filter Included
- Installation Hardware Included
- Add an Air/Water Separator \$199
- Limited 5 Year Warranty



## \$3595 2HP PRECISION VACUUM SYSTEM

- 2 HP Dual Solid State Motor
- 240Volt
- Machined Naval Brass Housing
- Spin Balance Impeller for Longevity
- Particulate Filter Included
- Installation Hardware Included
- Add an Air/Water Separator \$289
- Limited 5 Year Warranty




## \$4695 NEW! 4HP PRECISION VACUUM SYSTEM

- 4 HP Dual Solid State Motor
- 240Volt
- Machined Naval Brass Housing
- Spin Balance Impeller for Longevity
- Particulate Filter Included
- Installation Hardware Included
- Add an Air/Water Separator \$289
- Limited 5 Year Warranty



1 HP VACUUM AIR  
WATER SEPARATOR  
\$199



2 HP VACUUM AIR  
WATER SEPARATOR  
\$289



LOW VOLTAGE  
CONTROL BOX  
\$295



LOW VOLTAGE  
CONTROL PANEL  
\$225



WATER SOLENOID  
FILTER  
\$995

# DRY VACUUM SYSTEMS



## ROADRUNNER HYBRID 2.5HP VARIABLE SPEED

- 2.5 HP Variable Speed
- Real Time Sensor Adjust Vacuum Usage
- Noise level 40dB
- 240 Volt / Amps - 6.4
- Vacuum in / HG - 7 - 9
- Stainless Steel Tank
- Tank Wash Included
- Reduced Electrical Costs
- Virtually No Water Usage
- Environmentally Friendly
- 5 Year Limited Warranty

\$5295

## ROADRUNNER DUAL HYBRID 5HP VARIABLE SPEED

- DUAL 2.5 HP Variable Speed Motors
- Real Time Sensor Adjust Vacuum Noise level 40dB
- 240 Volt / Amps - 12.8
- Vacuum in / HG - 7 - 9
- Stainless Steel Tank
- Tank Wash Included
- Reduced Electrical Costs
- Virtually No Water Usage
- Environmentally Friendly
- 5 Year Limited Warranty

\$8995



## ROADRUNNER HYBRID 3.5HP VARIABLE SPEED

- 3.5 HP Variable Speed
- Real Time Sensor Adjust Vacuum Usage
- Noise level 40dB
- 240 Volt / Amps - 6.4
- Vacuum in / HG - 7 - 9
- Stainless Steel Tank
- Tank Wash Included
- Reduced Electrical Costs
- Virtually No Water Usage
- Environmentally Friendly
- 5 Year Limited Warranty

\$6295

## ROADRUNNER DUAL HYBRID 7HP VARIABLE SPEED

- DUAL 3.5 HP Variable Speed Motors
- Real Time Sensor Adjust Vacuum Noise level 40dB
- 240 Volt / Amps - 12.8
- Vacuum in / HG - 7 - 9
- Stainless Steel Tank
- Tank Wash Included
- Reduced Electrical Costs
- Virtually No Water Usage
- Environmentally Friendly
- 5 Year Limited Warranty

\$10195



# OPERATOR LIGHTS



## INFINITY DENTAL LIGHT INSTALLS ON ANY DENTAL CHAIR

- Light Head with 3 Axis Movement
- LED Light Source
- Adjustable Intensity
- Color Temperature: 4800 kelvin
- Color Rendering Index 92
- Red Coloring Index 60
- Focal Shape Oval
- Counter Balance, Drift Free Arm
- Dimension Specifications  
Length Horizontal Arm - 34"  
Length Flex Arm - 26"  
Flex Arm Vertical - 31"w

**\$1295**



## INFINITY LED CEILING LIGHT MOUNTS ON 8' TO 10' CEILINGS

- Light Head with 3 Axis Movement
- LED Light Source
- Adjustable Intensity
- Color Temperature: 4800 kelvin
- Color Rendering Index 92
- Red Coloring Index 60
- Focal Shape Oval
- Counter Balance, Drift Free Arm
- Dimension Specifications  
Length Horizontal Arm - 34"  
Length Flex Arm - 26"  
Flex Arm Vertical - 31"

**\$1695**



## DIAMOND DENTAL LIGHT INSTALLS ON ANY DENTAL CHAIR

- Light Head with 2 Axis Movement
- Xenon Quartz Halogen Bulb
- Dual Intensity - High/Low Illumination
- Color Temperature: 4200 to 5000 kelvin
- Ergonomic Handles
- Simple Lever Activation
- Glass Reflector w/ Shadow Free Illumination
- Counter Balance, Drift Free Arm
- Dimension Specifications  
Length Horizontal Arm - 34"  
Length Flex Arm - 26"  
Flex Arm Vertical - 31"

**\$895**



## DC - WALL MOUNTED INTRAORAL X-RAY MACHINE

- Light Weight Tube Head
- Easy to Position
- Computer Controlled Exposures
- Zero Drift Arm
- Digital and Film Compatible
- Customizable Installation (Wall or Cabinet\*)
- \*May require optional features
- High Frequency DC Tube Head

**Call for Pricing**



## AC - WALL MOUNTED INTRAORAL X-RAY MACHINE

- Light Weight Tube Head
- Easy Position
- Zero Drift Arm
- Manual Controlled Exposures
- Digital and Film Compatible
- Customizable Installation (Wall or Cabinet\*)
- \*May require optional features

**Call for Pricing**



EQUIPMENT ACCESSORIES



Electric Foot Control  
\$325



Articulating Headrest  
\$375



Computer Monitor Bracket  
\$250



Self Contained Water  
System \$225



Infinity Piezo Scaler Built-in  
\$595



Fiber Optic System  
\$395



Unit Chair Control Module  
\$295



Chair Swivel  
\$350



Utility Tray  
\$150



Wet/Dry Foot Control  
\$125



Water System Quick  
Disconnect \$75



Fiber Optic System  
\$395

EQUIPMENT ACCESSORIES



# STERILIZERS & ULTRASONICS



## TUTTNAUER EZ 9 Plus STERILIZER W/ Printer

- Chamber Size - 9" x 19.8"
- Cycle Time 24 Minutes (Hot Start)
- Auto Drying Cycle
- Low Water Sensor
- Preprogramed Sterilization Cycles
- 120 volt - 12 amp
- Easy Drain for Water Reservoir
- Tray Capacity - 3 Full

\$4995

TUTTNAUER CSU1  
ULTRASONIC  
\$595  
1 Gallon



## TUTTNAUER EZ 10 Plus STERILIZER W/ Printer

- Chamber Size - 11" x 19.8"
- Cycle Time 14 Minutes (Hot Start)
- Auto Drying Cycle
- Low Water Sensor
- Pre-Programmed Sterilization Cycles
- 120 volt - 12 amp
- Easy Drain for Water Reservoir
- Tray Capacity - 5 Full

\$5395

TUTTNAUER CSU3  
ULTRASONIC  
\$1395  
3 Gallon



TUTTNAUER  
DISTILLER  
\$490  
1 Gallon







Dental Chairs | Delivery Systems | Operator Seating | Vacuum Systems | Air Compressors | X-Rays | Sterilizers

Dansereau has the most comprehensive product warranty programs available for dentists today.  
We offer a 10 Year Warranty on our Dental Chairs & Delivery Systems with complete features listed below:

**5 Years 100% / 6-10 Years Prorated Annually**

Dental Chair Electric Motors, Chair Base Frame and Chair Frame

**3 Years 100% / 4-10 Years Prorated Annually**

Delivery System Internal Tubing (against leaks), Master Valve Regulators, Handpiece Control Blocks, Light Posts, Automatic Handpiece Holders, Delivery System Water and Air Flow Control Valves, Stool Cylinders, Manual Chair Switches, Delivery System Flex Arm, Delivery System Telescoping Arm, Cuspidor Frame.

**1 Year 100%**

Diamond Dental Light, Chair Foot Control, Stool Ring, Stool Seat, Stool Casters, Circuit Board and Prepositioning Switches, Upholstery\* (Limited\*)

**6 Months 100%**

Base and Unit Covers parts, 3-way syringes, Delivery System Foot Control and Fiber Optics Systems.

**30 Days 100%**

Dental Light Bulbs, Fiber Optic Light Bulbs, X-ray View Box Light Bulbs

**Broken or User Damaged Items are Not Covered under Warranty**

Minimum Maintenance is required in order to guarantee full use of warranties. All dental delivery systems require clean dry air to allow the proper lifespan to the micro valves within the dental delivery system. If your air compressor does not have an adequate drying and filtration system the warranties of the micro valves, on/of toggle valves, relays, regulators, master valves, foot controls, 3 ways syringes and automatic blocks may be void. \*Upholstery warranty does not include discoloration and damage as a result of improper use and maintenance.

**GUARANTEE**

Dansereau will guarantee against defects in material and workmanship. Defective parts will be repaired or replaced: no claim for labor or consequential damages will be allowed. All replacement parts claimed defective will be filled in the usual manner, such invoice being subject to adjustment after the alleged defective material has been returned and inspect by the manufacturer. The warranty does not cover damage resulting from improper installation, maintenance, accident or misuse.

**Terms: Prepaid/Financing is Available**

Freight: All Shipments are FOB Corona, CA unless otherwise specified by purchaser on the order. Dansereau is not responsible for any freight unless otherwise specified in writing.  
Returned Product: All orders are non-cancelable. No parts can be returned without a factory authorization. All parts are subject to 25% re-stocking fee and must be returned Prepaid Shipping.

(800) 423-5657  
[www.dhpdental.com](http://www.dhpdental.com)