




UNIVERSAL CART SYSTEM SPECS

Universal Cart Specs

- Height Adjustable U Cart 35" to 27"
- Rotating Unit Head move from Right to Left Side of Cart
- Rotating Assistants Package mover from Right to Left Side of Cart
- Complete Junction Box - Not Shown
- Laminated Worksurface 21" x 18" - can be altered by request.



NOTE: Verify if Self Contained Water System is being used as Public Water is not required

- Critical Notes:**
- 1) Verify Water Requirements as Dental Delivery Systems with Self Contained Water Systems will not require Water in J Box
 - 2) Verify Vacuum Model as Line Size May Change
 - 3) Uniform Plumbing Code requires any Dental Delivery System to installed to the Public Water System an Independent Backflow Device Verify with the Local Building and Safety

Finished Floor